18/1978

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STANDARD CONDITIONS FOR CALCIUM

STOCK STANDARD SOLUTION

Calcium, 500 µg/ml. To 1.249 g of primary standard calcium carbonate, CaCO₃, add 50 ml of deionized water. Add dropwise a minimum volume of HCl (approx. 10 ml) to effect complete solution of the CaCO₃. Dilute to 1 liter with deionized water.

OPERATING PARAMETERS

Perkin-Elmer Instrument Settings

Instrument	503, 403, 306, 305(A,B), 303	300, 300S	103, 107, 290(B)
Wavelength - 422.7 nm	211 - VIS	423	459
Slit Setting	4 (1.4 nm)	0.7 nm	0.7 nm
Light Source	Hollow Cathode Lamp		
Flame Type	Air-acetylene flame Reducing		

SENSITIVITY

For the standard conditions described above, the sensitivity is about 0.08 µg/ml Ca for 1% absorption. A standard containing 4 µg/ml Ca will typically give an absorbance reading of about 0.22 absorbance units (about 40% absorption).

LINEAR WORKING RANGE

For the standard conditions described above, the working range for Ca is linear up to concentrations of approximately 7 µg/ml in aqueous solution.

March 1973

Ca

Standard Conditions

-2-

OTHER ANALYTICAL LINES

<u>Wave-length</u>	503, 403, 306, 305(A,B), 303	300, 300S	103, 107, 290(B)	<u>Slit</u>	<u>Relative Sensitivity</u>
422.7 nm	211 - VIS	423	459	0.7 nm	1.
239.9 nm	240 - UV	240	129	0.7 nm	120.

FLAME ADJUSTMENT

The absorption of calcium is dependent on the fuel-to-air ratio and the height of the light beam above the burner. The fuel flow should be adjusted to yield maximum absorption for a calcium standard before beginning the analysis.

INTERFERENCES

Silicon, aluminum, phosphate, and sulfate depress the sensitivity for calcium. Sidiropoulos (1471) has reported similar depression of calcium sensitivity by titanium and zirconium. Lanthanum or strontium at concentrations of 0.1 to 1% can be added to sample and standards to control these interferences, and simultaneously, a slight ionization interference for calcium in the air-acetylene flame. Acid concentration of the standards should be matched with that of the samples.

OTHER FLAMES

The nitrous oxide-acetylene flame will provide a sensitivity of 0.05 µg/ml for 1% absorption. Ionization interferences should be controlled by adding a large amount of alkali salt to the samples and standards. The analysis appears to be free from chemical suppressions in the nitrous oxide-acetylene flame. It is probably preferable to determine calcium in the nitrous oxide-acetylene flame especially in samples containing large amounts of silica.

FLAME EMISSION

The most sensitive emission wavelength for calcium is at 422.7 nm. A nitrous oxide-acetylene flame is recommended. Calcium can also be determined at the 239.9 nm wavelength, but with reduced sensitivity. An air-acetylene flame may also be used with reduced sensitivity.

March 1973

Mg
Standard Conditions

STANDARD CONDITIONS FOR MAGNESIUM

STOCK STANDARD SOLUTION

Magnesium, 1000 µg/ml. Dissolve 1.000 g of magnesium ribbon in a minimum volume of (1+1) HCl. Dilute to 1 liter with 1% (v/v) HCl.

OPERATING PARAMETERS

Perkin-Elmer Instrument Settings

Instrument	503, 403, 306, 305(A,B), 303	300, 300S	103, 107, 290(B)
Wavelength-285.2 nm	285 - UV	285	210
Slit Setting	4 (0.7 nm)	0.7 nm	0.7 nm
Light Source	Hollow Cathode Lamp		
Flame Type	Air-acetylene flame Oxidizing (lean, blue)		

SENSITIVITY

For the standard conditions described above, the sensitivity is about 0.007 µg/ml Mg for 1% absorption. A standard containing 0.3 µg/ml Mg will typically give an absorbance reading of about 0.19 absorbance units (about 35% absorption).

LINEAR WORKING RANGE

For the standard conditions described above, the working range for Mg is linear up to concentrations of approximately 0.5 µg/ml in aqueous solution.

March 1973

Mg

Standard Conditions

-2-

OTHER ANALYTICAL LINES

<u>Wave-length</u>	<u>503, 403, 306, 305(A,B), 303</u>	<u>300, 300S</u>	<u>103, 107, 290(B)</u>	<u>Slit</u>	<u>Relative Sensitivity</u>
285.2 nm	285 - UV	285	210	0.7 nm	1.0
202.6 nm	203 - UV	203	064	0.7 nm	24.

INTERFERENCES

Silicon and aluminum depress magnesium absorption in the air-acetylene flame. The addition of lanthanum (0.1 to 1.0%) will generally remove these interferences, as will the use of a nitrous oxide-acetylene flame.

FLAME EMISSION

The most sensitive emission wavelength for magnesium is at 285.2 nm. A nitrous oxide-acetylene flame is recommended. Magnesium can also be determined at the 333.7 and 382.9 nm wavelengths, but with reduced sensitivity. An air-acetylene flame can also be used, with reduced sensitivity.